

# ViveST™: A REVOLUTIONARY AMBIENT TEMPERATURE STORAGE AND TRANSPORTATION DEVICE FOR UTILIZATION IN INFECTIOUS DISEASE TESTING

AMP 2013  
ID09

Anita McClernon, Andrew Freeman and Daniel McClernon  
bioMONTR® Labs, Research Triangle Park, NC

## Introduction

Infectious disease monitoring often requires collection sites to ship biological specimens to reference laboratories for molecular based testing. Standard specimen storage and transportation methods have several disadvantages including temperature control and low-volume requirements.

Studies discussed here-in evaluate the use of ViveST, a transformational dried ambient matrix, for storage and transport of infectious plasma samples used with downstream molecular based viral load assays.

## Methods

- For all testing described, infectious plasma (HCV, HIV-1 or HBV) was loaded onto ViveST, dried, capped and stored at ambient temperature.
- After recovery from ViveST, specimens were analyzed using the following commercially available FDA approved viral load assays:
  - RealTime HCV Assay (Abbott Molecular, Des Plaines, IL)
  - COBAS® TaqMan® HCV Test, v2.0 for use with the High Pure System (Roche Diagnostics, Indianapolis, IN)
  - RealTime HIV-1 Assay (Abbott Molecular, Des Plaines, IL)
  - COBAS® TaqMan® HIV-1 Test, v2.0 for use with the High Pure System (Roche Diagnostics, Indianapolis, IN)
  - RealTime HBV Assay (Abbott Molecular, Des Plaines, IL)
- To assess precision (inter- and intra-assay) for each analyte/assay combination, multiple aliquots of infectious plasma specimens (3 or 5 depending on assay/analyte combination) with varying viral loads (3 levels) were analyzed on 3 separate runs after storage (RT) on ViveST for up to 10 days.
- To assess linearity for each analyte/assay combination, a high titer infectious plasma sample was diluted in normal human plasma (7 levels). Each level was loaded onto ViveST in triplicate, stored for 7 days (RT), recovered and tested on a single run (n = 21 each analyte/assay combination). For comparison purposes, identical aliquots of frozen plasma were tested simultaneously.

## Results (HCV Viral Load)


### Abbott RealTime HCV

- Precision results are summarized in Table 1. HCV infectious samples processed through ViveST yielded reproducibility results with a standard deviation <0.10 LOG IU/mL (intra-assay) and <0.07 LOG IU/mL (inter-assay). The 95%CI were <±0.11 (intra-assay) and <±0.04 (inter-assay).
- Analysis of diluted samples from ~1 to ~6 LOG IU/mL demonstrated a direct proportional relationship between the dilution factor and the number of HCV copies reported with R<sup>2</sup>>0.99 (see Figure 1).

Table 1. ViveST\_Abbott HCV Intra-assay and Inter-assay Precision

Concentration	ViveST_Abbott HCV Intra-assay and Inter-assay Precision								
	Low			Mid			High		
Timepoint (Day)	1	2	3	1	2	3	1	2	3
Replicates (n)	3	3	3	3	3	3	3	3	3
Mean (LOG IU/mL)	2.21	2.18	2.22	3.71	3.63	3.63	5.16	5.03	5.10
Std Dev	0.08	0.03	0.10	0.04	0.05	0.02	0.01	0.05	0.02
95% CI	0.03	0.03	0.11	0.01	0.06	0.02	0.00	0.06	0.02

Figure 1. ViveST\_Abbott HCV Analytical Measurement Range




### COBAS® TaqMan® HCV Test, v2.0

- Precision results are summarized in Table 2. HCV infectious samples processed through ViveST yielded reproducible results with a standard deviation of <0.15 LOG IU/mL (intra-assay) and <0.11 LOG IU/mL (inter-assay). The 95%CI were <±0.17 (intra-assay) and <±0.07 (inter-assay).
- Analysis of diluted samples from ~4 to ~6 LOG IU/mL demonstrated a direct proportional relationship between the dilution factor and the number of HCV copies reported with R<sup>2</sup>>0.99 (see Figure 2).

Table 2. ViveST\_Roche HCV Intra-assay and Inter-assay Precision

Concentration	ViveST_Roche HCV Intra-assay and Inter-assay Precision								
	Low			Mid			High		
Run #	1	2	3	1	2	3	1	2	3
Days Stored	1	3	7	1	3	7	1	3	7
Replicates	3	3	3	3	3	3	3	3	3
Mean (LOG IU/mL)	3.65	3.59	3.48	4.17	4.22	4.06	4.48	4.49	4.38
Std Dev	0.15	0.04	0.07	0.13	0.04	0.07	0.01	0.03	0.09
95% CI	0.17	0.04	0.08	0.15	0.04	0.08	0.01	0.03	0.10

Figure 2. ViveST\_Roche HCV Analytical Measurement Range



## Results (HIV-1 Viral Load)


### Abbott RealTime HIV-1

- Precision results are summarized in Table 3. HIV-1 infectious samples processed through ViveST yielded reproducibility results with a standard deviation <0.11 LOG c/mL (intra-assay) and <0.10 LOG c/mL (inter-assay). The 95%CI were <±0.10 (intra-assay) and <±0.05 (inter-assay).
- Analysis of diluted samples from ~1 to ~7 LOG c/mL demonstrated a direct proportional relationship between the dilution factor and the number of HIV-1 copies reported with R<sup>2</sup>>0.99 (see Figure 3).

Table 3. ViveST\_Abbott HIV-1 Intra-assay and Inter-assay Precision

Concentration	ViveST_Abbott HIV-1 Intra-assay and Inter-assay Precision								
	Low			Mid			High		
Run #	1	2	3	1	2	3	1	2	3
Days Stored	3	7	10	3	7	10	3	7	10
Replicates (n)	5	5	5	5	5	5	5	5	5
Mean (LOG c/mL)	2.71	2.58	2.66	3.64	3.58	3.60	4.72	4.61	4.62
Std Dev	0.08	0.08	0.11	0.06	0.03	0.06	0.03	0.09	0.05
95% CI	0.05	0.07	0.10	0.05	0.03	0.05	0.02	0.08	0.04

Figure 3. ViveST\_Abbott HIV-1 Analytical Measurement Range




### COBAS® TaqMan® HIV-1 Test, v2.0

- Precision results are summarized in Table 4. HCV infectious samples processed through ViveST yielded reproducible results with a standard deviation of <0.17 LOG c/mL (intra-assay) and <0.25 LOG c/mL (inter-assay). The 95%CI were <±0.19 (intra-assay) and <±0.16 (inter-assay).
- Analysis of diluted samples from ~1 to ~7 LOG c/mL demonstrated a direct proportional relationship between the dilution factor and the number of HCV copies reported with R<sup>2</sup>>0.99 (see Figure 4).

Table 4. ViveST\_Roche HIV-1 Intra-assay and Inter-assay Precision

Concentration	ViveST_Roche HIV-1 Intra-assay and Inter-assay Precision								
	Low			Mid			High		
Run #	1	2	3	1	2	3	1	2	3
Days Stored	1	4	7	1	4	7	1	4	7
Replicates	3	3	3	3	3	3	3	3	3
Mean (LOG c/mL)	2.49	1.97	2.06	3.16	2.86	2.90	5.04	4.61	4.72
Std Dev	0.01	0.04	0.12	0.07	0.06	0.09	0.02	0.17	0.05
95% CI	0.01	0.04	0.14	0.08	0.07	0.10	0.02	0.19	0.06

Figure 4. ViveST\_Roche HIV-1 Analytical Measurement Range



## Results (HBV Viral Load)


### Abbott RealTime HBV

- Precision results are summarized in Table 5. HBV infectious samples processed through ViveST yielded reproducibility results with a standard deviation <0.13 LOG IU/mL (intra-assay) and <0.10 LOG IU/mL (inter-assay). The 95%CI were <±0.14 (intra-assay) and <±0.06 (inter-assay).
- Analysis of diluted samples from ~1 to ~7 LOG IU/mL demonstrated a direct proportional relationship between the dilution factor and the number of HBV copies reported with R<sup>2</sup>>0.99 (see Figure 5).

Table 5. ViveST\_Abbott HBV Intra-assay and Inter-assay Precision

Concentration	ViveST_Abbott RealTime HBV Intra-assay and Inter-assay Precision								
	Low			Mid			High		
Run #	1	2	3	1	2	3	1	2	3
Days Stored	1	3	7	1	3	7	1	3	7
Replicates	3	3	3	3	3	3	3	3	3
Mean (LOG IU/mL)	3.61	3.57	3.66	4.74	4.70	4.83	5.81	5.82	5.84
Std Dev	0.04	0.05	0.06	0.02	0.13	0.10	0.02	0.04	0.03
95% CI	0.04	0.05	0.07	0.02	0.14	0.11	0.03	0.04	0.03

Figure 5. ViveST\_Abbott HBV Analytical Measurement Range



## Conclusions

- These studies demonstrate ViveST's utility for storing and transporting infectious plasma for subsequent RNA and DNA viral target analysis using various commercially available molecular based assays.
- ViveST provides significant cost savings as compared to shipping samples on dry ice and can enhance access to healthcare globally while significantly reducing the burden associated with shipping frozen samples.
- ViveST provides superior reproducibility and accuracy for HCV, HBV and HIV-1 viral load testing.

## Acknowledgments

The authors would like to thank ViveBio, LLC for their support during this study.

Send correspondence to:  
Anita McClernon  
amccclernon@biomontr.com

